

## Direct ELISA Protocol

### **Buffers and reagents**

Antigen or antibody should be diluted in coating buffer to immobilize them to the wells: PBS pH 7.4. See Gibco recipe.

### **Blocking solution**

Commonly used blocking agents are 1% BSA, serum, non-fat dry milk, casein, gelatin in PBS. – 1% BSA in TBS-T (tris- buffered saline with 0.05% Tween 20)

### **Rinse Solution**

1X TBS → no tween

### **Wash solution**

Usually PBS or Tris-buffered saline (pH 7.4) with detergent such as 0.05% (vol/vol) Tween20 (TBST).

### **Antibody dilution buffer**

Primary and secondary antibody should be diluted in 1% BSA in TBS-tween blocking solution to reduce non-specific binding.

Antibody should be stored in 2-8 degrees C; see minifridge.

### **Stop Solution**

4M HCL

---

Fisher:

- 96 well plates - we use Thermo Scientific Cat# 12-565-136
- Coating buffer- Gibco PBS pH 7.4 - Cat# 0-010-031

Sigma:

- Tween 20 - Cat# P7949-100mL
- BSA- Cat# A7906

Biologend:

- Substrate solution- TMB substrate set Cat# 421101
- Alternative coating buffer - Sodium bicarbonate buffer pH 9.5 Cat# 421701

Biorad:

- Tris Buffered Saline- Cat# 1706435



- Remove Wash Buffer after the third wash and **rinse** twice with 200µL per well of 1X TBS (no tween). Aspirate or decant and tap the plate gently over a stack of tissues to remove excess liquid.

### Detection

- TMB substrate comes in Substrate A and Substrate B. Mix the two substrate solutions in a conical tube with a 1:1 ratio. Calculate 100 uL x # of wells (round up the nearest mL so you don't run out. THIS IS CRUCIAL.
- Dispense 100 µl of the **TMB substrate solution** per well with a multichannel pipet.
  - o BE SURE TO CHANGE TIPS
  - o MAKE SURE NO BUBBLES IN THE TIPS IN THIS STEP
  - o ALSO BE FAST, AS COLOR WILL DEVELOP
- After substrate has been added to wells, place plate in a drawer so the color can develop in the dark
- After sufficient color development (if it is necessary) add 100 µl of stop solution to the wells
  - o APPLY TO WELLS IN THE SAME ORDER THAT SUBSTRATE WAS APPLIED.
  - o MUST CHANGE TIPS HERE.
- Read the absorbance (optical density) of each well with the Tecan Plate Reader.
  - o 450 nm filter absorbance.
  - o There is a file in the BIOMOD folder called ELISA.
  - o Contingency if you can't find the folder: The only change is the type of plate we use and the filter absorbance.
- Note: some enzyme substrates are considered hazardous (potential carcinogens), therefore always handle with care and wear gloves. ☺